AP16 IF / AP16 IF BLOT ELITE

Fast, Precise, Accurate

16 IFA slides Automated Processor

IF and BLOT fully automation

16 IFA slides and 24 BLOT strips

AP16 IF ELITE is an innovative and fast processor able to run 16 IFA slides with any kind of thickness, mask, well shape and number of rows. It is equipped with a sliding-sample barcode reader and a camera for slide barcode reading. Reliable and manufactured with the latest technologies, it is able to carry out: sample dilution, sample and reagent dispense, different incubation times and slide well different washing types with accuracy and precision.

AP16 IF BLOT ELITE is a powerful combo system for IFA and BLOT methods processing thanks to the easy and fast work area configuration change. In addition to all the AP16 IF ELITE features a dedicated Blot Module and camera allows to carry out a full automation for Dot/Western Blot, starting from sample dilution, sample and reagent dispense, incubation, strip washing up to the dot imaging processing and assessment for each strip. The software allows to edit, store and record patient data and reports.

APIG IF / APIG IF BLOT **ELITE**

Built-in barcode reader for samples, controls and reagents; BCR for slides (optional)

Average throughput 192 wells in a run with incubation times of 30 + 30 minutes

Sample trays 180 sample tubes Ø13 mm or 135 sample tubes Ø 16 mm; combination of them can be used in the same run

Supplied trays 264 positions for dilution cups; 16 positions for IFA slides; 20 controls; 8 conjugates (7 in case of Blue Evans dye) and 1 cleaning solution for the probe; 2 diluents for samples

External liquid tanks 2 washing buffers for slides, Wash and Waste tanks, all equipped with level sensors

Slide washing wells are washed one by one to avoid carryover

Washing modes multi-well, continuous flow or drop by drop for specific washing effectiveness

Liquid handling 2 probes: one to dispense and dilute, the second to aspirate from wells

Sample dilution time 20s with 2ml of washing solution

Sample dispense time 11s with 2 ml of washing solution

Dispense CV <1%, inaccuracy 1% at 100 µl

Aspiration CV <3%, inaccuracy 1% at 100 µl; minimum volume 5 µl, 1 µl resolution

Syringes 1000 µl volume, 1 µl resolution

Operating Functions Management of 8 methods in the same run. Definition of incubation times and volume for sample and reagents dispense. Different incubation times for each method. Easy and fast settings for new protocols.

Connection LIS interface capability, work list import from Host Computer; a single PC can run more instruments at the same time

PC specifications Intel I5 processor; 4 GB RAM; VGA 1024 x 768; HD with 100 GB free space; CD-ROM reader; USB 2.0 port; RS232 Serial Port

Operating system Windows 7 or Windows10

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