
# AP32 IF ELITE

Fast, Precise, Accurate  
40 IFA slides Processor

**AP32 IF ELITE** is an innovative system able to process 40 IFA slides. A special three needle system expedites the work session. The system uses a powerful proprietary software that allows to run more DAS' instruments at the same time.

It processes simultaneously slides with any kind of thickness, mask, well shape and number of rows, keeping the well known accuracy of all IF processors of DAS. The instrument is equipped with sample barcode reader; barcode reader for slides is available as well.

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# AP32 IF ELITE

<b>Built-in barcode reader</b>	for samples and for slides
<b>Average throughput</b>	320 wells in a run with incubation times of 30 + 30 minutes
<b>Sample trays</b>	250 sample tubes Ø12 mm or 190 sample tubes Ø 16 mm; combination of them can be used in the same run
<b>Supplied trays</b>	528 positions for dilution cups; 32 positions for IFA slides (optional 40 positions); 25 controls; 11 conjugates (10 in case of Blue Evans dye) and 1 cleaning solution for the probe; 6 diluents for samples
<b>External liquid tanks</b>	2 washing buffers for slides, Wash and Waste tanks, all equipped with level sensors
<b>Slide washing</b>	wells are washed one by one to avoid carryover
<b>Washing modes</b>	multi-well, continuous flow or drop by drop for specific washing effectiveness
<b>Liquid handling</b>	3 probes: two to dispense and dilute, the third to aspirate from wells
<b>Sample dilution time</b>	15s with 2ml of washing solution
<b>Sample dispense time</b>	6s (average) with 2 ml of washing solution
<b>Dispense</b>	CV <1%, inaccuracy 1% at 100 µl
<b>Aspiration</b>	CV <3%, inaccuracy 1% at 100 µl; minimum volume 5 µl, 1 µl resolution
<b>Syringes</b>	2 with 1000 µl volume each, 1 µl resolution
<b>Operating Functions</b>	Management of 10 methods in the same run. Definition of incubation times and volume for sample and reagents dispense. Different incubation times for each method. Easy and fast settings for new protocols.
<b>Connection</b>	LIS interface capability, work list import from Host Computer; a single PC can run more instruments at the same time
<b>PC specifications</b>	Intel I5 processor; 4 GB RAM; VGA 1024 x 768; HD with 100 GB free space; CD-ROM reader; USB 2.0 port; RS232 Serial Port
<b>Operating system</b>	Windows 7 or Windows10
<b>Power requirements</b>	230 or 115 V AC, 50 - 60 Hz, 450VA
<b>Measurements</b>	113 x 77 x h 75 cm
<b>Weight</b>	110 Kg.

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The shown features could be altered without notice

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